

Hawk-i 5.1.0.6 Release Notes

Hawk-i is DICOM viewer and utility software designed for viewing multi-modality medical images especially for cardiovascular multi-frame movies.

What's New

A few of the features new to this release include:

- CD-ROM button (Play DICOM CD) Improvements

Enable /Disable this button if CD-ROM drive is power ON/OFF while Hawk-i is running.

- Sidebar docking Improvements
- Bug fixes

1)Failed to view 16bit MR image which column is 1022.

2)Failed to save as movie from the DICOMCD which comes from some particular vendors.

3)Failed to adjust CT Level/Width.

- More stable and faster

With every release, Hawkie works to make Hawk-i more stable and faster.

- Other Bug Fixes

For brevity some small bugs are not reported :-)

Install and Uninstall

Once you have gotten the Hawk-i newer setup binary, just run the setup simply. If you prefer to install old version of Hawk-i, please uninstall it beforehand .

To uninstall HAWk-i, simply run Start-All Program-Hawkie-Uninstall

Feedback and Bug Reports

If you find a problem that you think might be a bug in Hawk-i, please contact us.

We want to hear your feedback about Hawk-i. Please feel free to contact us .

Hawkie MediTech, Inc.

Hawkie Makes IT Simple...

Email: wu@hawkie.net url: <http://hawkie.net>